

## *SURFACE ROUGHNESS AND FEATURES*

— OF A —

## *SOLAR CELL*



Prepared by  
**Xavier Herrera-Keehn, Jocelyn Esparza, and Pierre Leroux**



## Introduction


Maximizing a solar cell's energy absorption is key for the technology's survival as a renewable resource. The multiple layers of coating and glass protection allow for the absorption, transmittance, and reflection of light that is necessary for the photovoltaic cells to function. Given that most consumer solar cells operate at 15-18% efficiency, optimizing their energy output is an ongoing battle.

## Importance of Solar Panel Testing

Studies have shown that surface roughness plays a pivotal role in the reflectance of light. The initial layer of glass must be as smooth as possible to mitigate reflectance of light, but the subsequent layers do not follow this guideline. A degree of roughness is necessary at each coatings interface to another to increase the possibility of light scattering within their respective depletion zones and increase the absorption of light within the cell<sup>1</sup>. Optimizing the surface roughness in these regions allows the solar cell to operate to the best of its ability and with the Nanovea HS2000 High Speed Sensor, measuring surface roughness can be done quickly and accurately.

# Measurement Objectives

In this study we will display the capabilities of the Nanovea HS2000 High Speed Sensor by measuring the surface roughness and geometric features of a photovoltaic cell. For this demonstration a monocrystalline solar cell with no glass protection will be measured but the methodology can be used for various other applications.



**Solar Cell sample on Nanovea HS2000L Profilometer**

# Measurement Parameters

The following test parameters were used to measure the surface of the solar cell.

Solar Cell	
Optical Pen	L1
Acquisition rate	100 pps
Averaging	1
Measured surface	30mm x 10mm
Step size	5 $\mu$ m x 5 $\mu$ m
Measurement Time (h:m:s)	00:13:43


Table 1: Test parameters used




Sample of solar cell analyzed

# Results and Discussion

Depicted below is the 2D false color view of the solar cell and an area extraction of the surface with its respective height parameters. A gaussian filter was applied to both surfaces and a more aggressive index was used to flatten the extracted area. This excludes form (or waviness) larger than the cut-off index, leaving behind features that represent the solar cell's roughness.



**Figure 1: False color view of the photovoltaic cell**



**Figure 2: False color view of the photovoltaic cell**

# Results and Discussion




Figure 3: Surface roughness of the first layer of the photovoltaic surface




Figure 4: 3D view of the first layer of the photovoltaic surface

# Results and Discussion




Figure 5: Profile roughness amplitude parameters of the first layer of the photovoltaic surface

A profile was taken perpendicular to the orientation of the gridlines to measure their geometric characteristics which is shown below. The gridline width, step height, and pitch can be measured for any specific location on the solar cell.




Figure 6: Width measurement of the gridline

# Results and Discussion





Figure 7: Pitch measurement of the gridline



Parameters	Unit	Step 1	Step 2	Step 3
Maximum height	µm	10.48	10.24	12.30
Mean height	µm	9.130	8.972	11.75

Figure 8: Step height measurement of the gridline



## Conclusion


In this study we were able to display the Nanovea HS2000 Line Sensor's ability to measure a monocrystalline photovoltaic cell's surface roughness and features. With the ability to automate accurate measurements of multiple samples and set pass fail limits, the Nanovea HS2000 Line Sensor is a perfect choice for quality control inspections.

To learn more about [Nanovea Profilometer](#) or [Lab Services](#).

## Reference

- 1 Scholtz, Lubomir. Ladanyi, Libor. Mullerova, Jarmila. "Influence of Surface Roughness on Optical Characteristics of Multilayer Solar Cells" *Advances in Electrical and Electronic Engineering*, vol. 12, no. 6, 2014, pp. 631-638.

## HS2000 Profiler



### High Speed Inspection & Precision Flatness Measurement

Advanced Automation with customizable options


High Speeds up to 384,000 points per second

Designed for large area flatness measurement

Full granite base with integrated anti-vibration table

[Learn More about the HS2000](#)

## Jr25 Profiler



### High Quality with Portability and Versatility

Full portability with Nanometer Resolution

25mm XY Scan Area

White-light Chromatic Confocal Technology

No Sample Size or Weight Limits

[Learn More about the Jr25](#)

## Thank you for reading!

We appreciate your interest in our technology and services. Read more about all of our product line and lab services at [www.nanovea.com](http://www.nanovea.com)

Call to Schedule a demo today!

If you have any questions please email us at [info@nanovea.com](mailto:info@nanovea.com)

## Recommended Reading

Check out our other application note :

[Nano Mechanical Characterization  
of Spring Constants](#)

**NANOVEA**  
A Better Measure

www.nanovea.com  
info@nanovea.com  
euroinfo@nanovea.com  
mexinfo@nanovea.com  
(949) 461-9292

**NANO MECHANICAL CHARACTERIZATION**  
— OF —  
**SPRING CONSTANTS**



Prepared by  
**Xavier Herrera-Keehn & Jocelyn Esparza**

Today's standard for tomorrow's materials. © JULY 2019 NANOVEA