

**Wafer Coating Thickness Measurement  
Using 3D Profilometry**



Prepared by  
**Jocelyn Esparza**

6 Morgan, Ste156, Irvine CA 92618 · P: 949.461.9292 · F: 949.461.9232 · [nanovea.com](http://nanovea.com)


**Today's standard for tomorrow's materials. © 2018 NANOVEA**

## INTRO:

Silicon wafers are widely used in the making of integrated circuits and other micro devices used in a vast number of industries. A constant demand for thinner and smoother wafers and wafer coatings makes the Nanovea 3D non-contact profiler a great tool to quantify coating thickness and roughness of just about any surface. The measurements in this article were taken from a coated wafer sample in order to demonstrate the capabilities of our 3D Non-Contact Profilometer.

### IMPORTANCE OF 3D NON CONTACT PROFILOMETER FOR STEP HEIGHT AND ROUGHNESS MEASUREMENT

Unlike other techniques such as touch probes or interferometry, the 3D Non-Contact Profilometer, using axial chromatism, can measure nearly any surface, sample sizes can vary widely due to open staging and there is no sample preparation needed. Nano through macro range is obtained during surface profile measurement with zero influence from sample reflectivity or absorption, has advanced ability to measure high surface angles and there is no software manipulation of results. Easily measure any material: transparent, opaque, specular, diffusive, polished, rough etc. The technique of the Non-Contact Profilometer provides an ideal, broad and user friendly capability to maximize surface studies when step height will be needed; along with the benefits of combined 2D & 3D capability.




### MEASUREMENT OBJECTIVE

In this application the nanovea PS50 is used to measure the step-height of a silicon wafer coating and its roughness. Assuming that the coating thickness is relatively uniform, the coating was removed to expose the substrate in order to measure it as the reference plane. The height of the coating will be measured to characterize the thickness and roughness of the coating.


## RESULTS: 3D Surface

The 3D View and False Color View of the selected area on the sample surface after the form has been removed. It provides users a straightforward tool to directly observe the distribution and morphology of the surface features from different angles.

3D View




False Color



## RESULTS: 2D Surface Analysis



The surface features can be quantified in numerous ways using built-in software tools as shown in the following figures as examples. It can be observed that the average coating thickness is 3.532  $\mu\text{m}$  and has an average surface roughness of 0.111  $\mu\text{m}$

Area  
Step Height



Differential parameters	P2 - P1	Unit
Zmean(higher) - Zmean(lower)	3.485	$\mu\text{m}$


## Profile Step Height



Parameters	Unit	Step 1	Step 2
Maximum height	µm	3.942	4.122
Mean height	µm	3.543	3.340

## Surface Area Roughness

ISO 25178	
Height Parameters	
Sa	0.1113 µm
Sz	2.690 µm
Sv	1.515 µm
Sp	1.175 µm
Sq	0.1417 µm
Ssk	0.1406
Sku	4.189



## CONCLUSION:

In this application, we have shown how the Nanovea 3D Non Contact Profilometer is capable of precisely characterizing the height and surface roughness of a metal coating. The data shows the coating to have a mean step-height of 3.532 µm and average surface roughness of 0.111 µm. The data shown here represents only a portion of the calculations available in the analysis software.

Learn more about the [Nanovea Profilometer](#) or [Lab Services](#)