




NANOVEA

OPTICAL PROFILER



Offering more than **25 Years** of Material Science Experience



RESEARCH AND CONSULTATION

Extensive range of research content such as brochures, application notes, publications, and videos.



EXPERT ASSISTANCE

Dedicated Profilometry experts happy to guide you through any question or project request.



CUTTING EDGE INNOVATION

At Nanovea we are always developing cutting edge technologies and standards. We innovate our instruments so that you can innovate your own products.



PRE AND POST INSTALLATION SUPPORT

Full walk-through and guide to make sure the instrument is installed perfectly. Dedicated support team to help you after your instrument has been installed.

INSTRUMENTS



ST400 OPTICAL PROFILER

- 200 x 150mm XY stages
- Video imaging integration
- Ideal for wide range of samples with varied geometries
- Chromatic confocal sensors w/ speed up to 200 times faster
- Rotational stage parallel or perpendicular to the testing plate
- Height sample clearance up to 200mm



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

THE STANDARD FOR PROFILOMETRY



X - Y SCAN AREA
200 x 150mm Motorized



HEIGHT RANGE
2.5nm to 25mm



DESKTOP DIMENSIONS
62 x 62 x 82cm



SCAN SPEED
40mm/s

ST500 LARGE AREA OPTICAL PROFILER

- High speed large area measurement w/ high speed sensor
- 400 mm XY axis travel with a maximum speed up to 200 mm/s
- Video zoom camera to provide automated functions
- Measurements with a user friendly desktop platform



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

HIGH SPEED AND LARGE AREA MEASUREMENT



X - Y SCAN AREA
400 x 400mm Motorized



HEIGHT RANGE
2.5nm to 25mm



DESKTOP DIMENSIONS
97 x 72 x 92cm



SCAN SPEED
200mm/s

JR25 PORTABLE OPTICAL PROFILER

- First truly portable non contact profilometer
- Weight less than 5.5 kg
- Lab quality results on the go
- Measurement capabilities up to 25mm x 25mm
- Able to measure samples at difficult angles
- Possible integration into automated robot arms and other equipment



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS




Atomic Force Microscope




Zoom Microscope

LABORATORY QUALITY RESULTS IN ANY LOCATION



X - Y SCAN AREA
25 x 25mm Motorized



HEIGHT RANGE
2.5mm to 25mm




DESKTOP DIMENSIONS
20 x 30 x 17cm



SCAN SPEED
20mm/s

PS50 COMPACT OPTICAL PROFILER

- Most advanced compact profilometer
- Small and simple footprint
- Measurement capabilities up to 50mm x 50mm
- All testing capabilities in compact version



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

MOST ADVANCED **COMPACT BENCHTOP**



X - Y SCAN AREA
50 x 50mm Motorized



HEIGHT RANGE
2.5nm to 25mm



DESKTOP DIMENSIONS
38 x 33 x 43cm



SCAN SPEED
20mm/s

JR100 PORTABLE & HIGH SPEED OPTICAL PROFILER

- Fast measurement (without stitching) using a 100 mm XY axis travel
- Z stage allows setup of measurements at various starting heights
- A high speed sensor gives ultra fast measurements at 382,000 points per second.
- Powerful for quality control



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

PORTABILITY AND HIGH SPEED



X - Y SCAN AREA
100 x 100mm Motorized



HEIGHT RANGE
2.5mm to 25mm



DESKTOP DIMENSIONS
44 x 49 x 32cm



SCAN SPEED
20mm/s

AFMPRO OPTICAL PROFILER

- 150 x 200mm XY stages and an adjustable height clearance of up to 140mm
- High magnification microscopy
- AFM expands the 3D capabilities into the sub nanometer range
- AFM gives the best lateral accuracy compared to optical techniques
- Easy to select zones on the video to be scanned



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

OPTICAL PROFILER *WITH AFM MODULE*



X - Y SCAN AREA
200 x 150mm Motorized



HEIGHT RANGE
2.5nm to 25mm




DESKTOP DIMENSIONS
64 x 64 x 82cm



SCAN SPEED
40mm/s

HS2000 BEST STABILITY OPTICAL PROFILER

- Granite base provides superior stability
- Automated inspection for quality control
- Workstation included to create fully contained stand alone instrument
- Excellent for roughness measurements, combined with advanced automation features



3D CHROMATIC SENSORS



Standard Sensor



High Speed Sensor

VIDEO OPTIONS





Atomic Force Microscope




Zoom Microscope

HIGH SPEED AND PRECISION FLATNESS TOOL



X - Y SCAN AREA
400 x 500mm Motorized



HEIGHT RANGE
2.5nm to 25mm



DESKTOP DIMENSIONS
101 x 106 x 195cm



SCAN SPEED
500mm/s


TECHNIQUE



CHROMATIC CONFOCAL

Chromatic Confocal technique uses white light that passes through a series of lenses with high degree of chromatic aberrations. Each wavelength will focus at a different distance creating the vertical measurement range. When a surface of interest is within the measurement range a single wavelength of the white light will be in focus while all others will be out of focus.

Only the focused wavelength will pass through the pin hole filter to reach the CCD spectrometer. The physical wavelength measured corresponds to a vertical position.



◆ NO USE OF COMPLEX ALGORITHMS

◆ NO LEVELING REQUIRED

LATERAL RESOLUTION vs ACCURACY

THE PROBLEM WITH OTHER TECHNIQUES



THEM

Camera pixel size or display resolution size are often used as lateral resolution to impress clients. For these, complex algorithms used to determine what is actually in focus gives a very different story of actual accuracy especially on complex surfaces.


US

Chromatic Confocal lateral accuracy is determined by physics and directly related to the spot size of the light.

LASER SCANNING CONFOCAL MICROSCOPE

VS

WHITE LIGHT CHROMATIC CONFOCAL



LASER RADIATION

Laser Light Health Hazard
Need for care of reflected light

Change in wavelength of laser light
affects results on the same sample

Non significant “display resolution”
Height & lateral accuracy fixed by objective used
Complex accuracy calculations

Alpha blending algorithms to combine
layer by layer data for complex accuracy calculation

Limited fixed field of view
Inaccurate stitching algorithms for larger surfaces

Data Acquisition speed 7900 Hz

SAFE WHITE LIGHT

UNIFORM BROAD WHITE LIGHT SPECTRUM
No effect of light intensity on results

INDEPENDENT LATERAL & HEIGHT ACCURACY
Any scan area at selected height accuracy

NO ALGORITHMS
Physical Wavelength Measured = Accurate Height

NO STITCHING
Continuous scanning of larger surfaces
Accuracy constant across any measurement size

50x FASTER
High Speed Sensor 384000 Hz

SCANNING A COIN

50x OBJECTIVE vs HIGH SPEED SENSOR (950 μm)

LATERAL ACCURACY

For 50x objective (370 x 277 μm)

$\pm 2\%$ of measuring value

$\pm 2\% \times 370 \mu\text{m}$

$\approx 15 \mu\text{m}$

w/ stitching algorithms $\gg 15 \mu\text{m}$



Step size:

$= 5 \mu\text{m}$

ULTIMATE LIMIT: $0.9 \mu\text{m}$

3x BETTER LATERAL ACCURACY

HEIGHT ACCURACY

$\approx 0.2 + L/100 \mu\text{m}$

$\approx 0.2 + 950/100 \mu\text{m}$

$\approx 9.7 \mu\text{m}$



950 μm range

$\approx 0.6 \mu\text{m}$

ULTIMATE LIMIT: $0.014 \mu\text{m}$

16x BETTER HEIGHT ACCURACY

AREA TESTED

STITCHING REQUIRED

scans (25 x 25 mm)

$25\,000 \mu\text{m} / 370 \mu\text{m} \times 25\,000 \mu\text{m} / 277 \mu\text{m}$

68×91

$= 6188 \text{ scans}$



NO STITCHING

Constant accuracy across any measurement size

1 SCAN


TEST TIME

6 sec per scan

+ 4 sec displacement & stitching

$= 10 \text{ sec/scan} \times 6188 \text{ scans}$

$= 61860 \text{ seconds} (\approx 17 \text{ hours})$




Scan time (25 x 25 mm)

$= 29.6 \text{ seconds}$


2090x FASTER

TESTING SOLUTIONS



ROUGHNESS | FINISH

- One second Ra measurement
- Any materials or surface complexity (3D or 2D)
- Automotive roughness finish standards





ROUGHNESS | FINISH ANALYSIS

- Ra | Sa profile & surface average roughness
- Rq | Sq profile & surface rms roughness
- Rz | Sz maximum height
- Sp | Sv maximum peak & pit height
- SKu | Ssk kurtosis & skewness of height distribution
- Bearing ratio and index
- Sk kernel roughness depth
- Spk | Svk reduced peak height & valley depth
- Sr1 | Sr2 upper & lower material ratio
- Sci & Svi core & valley fluid retention index


AND MORE

TEXTURE

- Isotropic & anisotropic surfaces
- Hills and valleys analysis



Texture Direction




TEXTURE ANALYSIS

- % of isotropy
- 1st, 2nd and 3rd direction
- % of periodicity
- Period
- Density of peaks
- Peak curvature (pointed or rounded)
- Average area of valleys & hills
- Average volume of valleys & hills

AND MORE

FLATNESS | WARPAGE

- Flatness <1µm over 500mm with no correction




FLATNESS | WARPAGE ANALYSIS

- 3D & 2D surface waviness & flatness
- FLTt peak to valley flatness deviation of the surface
- Best polynomial match
- FLTp peak to reference flatness deviation
- Material & bearing ratios
- FLTv reference to valley flatness deviation
- Distance measurement
- FLTq rms flatness deviation

AND MORE

VOLUME | AREA

- Surface subtraction & volume lost
- Corrosion analysis
- Motif and grain analysis





VOLUME | AREA ANALYSIS

- Volume of void, hills or valleys
- # of grains & average size
- Sdar | Spar developed surface area & projected area
- Area & perimeter of grains
- Volume of void & material from given height
- Height, area, volume of motifs
- Map area above or below given heights (%µm²)
- Max and min pitch of motifs
- Mean thickness of void & material from given height

AND MORE

GEOMETRY AND SHAPE

- Direct comparison to CAD geometry
- Curvature, radius, angles
- Lateral dimension
- Drill bit studies
- Cutting tools studies





GEOMETRY AND SHAPE ANALYSIS

- | | |
|------------------------------|--|
| • Radius of curvature | • Contour analysis |
| • Relative angle measurement | • Rake and wedge angle of drill bit |
| • Distance measurement | • K symmetry of cutting edge |
| • Mean diameter | • S alpha and gamma dist apex to end of clearance & rake roundness |

AND MORE

STEP HEIGHT | THICKNESS

- Measure through transparent materials
- Transparent film and coating thickness down to 20nm
- Steps from 20nm to 25mm



STEP HEIGHT | THICKNESS ANALYSIS


- | | |
|-------------------------------------|--------------------------------|
| • Point to point | • 3D or 2D map of thickness |
| • Point to plane | • Thickness distribution curve |
| • Maximum, minimum and mean heights | |

AND MORE


MICROSCOPE VIDEO IMAGING

Available on : **ST400, ST500, AFMPRO, & HS2000**

- Ultra zoom lens with coax lighting & detent
- Large area stitching capability
- Color video camera (1200x1600)
- Maximum magnification of 8000X
- Three positions turret (optional)



Broadview map selection tool



PRVision for machine vision capability






Image area selection measurement and image overlay



ATOMIC FORCE MICROSCOPE

Available on : **AFMPRO**

- Scan of XY 110 μm | high resolution XY 25 μm
- Lateral resolution 1.7nm
- Static, dynamic and extended modes
- Max Z range 22 μm | 5 μm
- Integrated video camera
- AFM to/from indenter position or video imaging with accuracy of < 0.2 μm



ADVANCED AUTOMATION



- Automatic focus (optical and microscope)
- Automatic analysis template
- Multi sample handling macros
- Easy selection of area under microscope
- Automatic dual frequency for surfaces with varying reflections
- Custom mounting setup of sensors for inline roughness QC
- Pattern recognition, database communications, pass/fail limits



VISIT OUR APPLICATION NOTES LIBRARY

nanovea.com/app-notes

Nanovea Optical Profilers **measure any material with a wider range of measurement** than any other Profilometer.



BASE	Jr25	Jr100	PS50	ST400	ST500	HS2000
Type	Portable	Portable & Fast	Compact	Standard	Large Area	Zero Noise / Flatness
X-Y Stage Travel	25 x 25mm	100 x 100mm	50 x 50mm	200 x 150mm	400 x 400mm	400 x 500mm
Z Axis	30mm Manual	25mm Manual	30mm Manual	50mm Motorized	50mm Motorized	100mm Motorized
Maximum X-Y Speed	20 mm/s	20 mm/s	20mm/s	40mm/s	200mm/s	500mm/s
System Dimensions	20 x 30 x 17cm	44 x 49 x 32cm	38 x 33 x 43cm	62 x 62 x 82cm	97 x 72 x 92cm	101 x 106 x 195cm
Rotational Options	N/A	N/A	N/A	Stage or Cylinder	Stage or Cylinder	Software
Video Microscope	N/A	N/A	N/A	Available	Available	Available
Max Sample Weight	No Limit	No Limit	8Kg	8Kg	4-8Kg	7Kg
High Speed Line Sensor	N/A	Included	N/A	Available	Available	Available
Customizable	50mm Stage Travel	N/A	N/A	4 axis & AFM	4 axis	5 axis

MEASUREMENT TECHNOLOGY

Technique — Non Contact • Chromatic Light
 Data Stitching — Not Required within X-Y Stage Travel
 Materials Types — ALL - Including Dark, Transparent, & Reflective
 Max Surface Angle — Up To 87°
 Max Vertical Resolution — 1nm

More Information at
nanovea.com/profilometers

STANDARD SENSOR (Single Point)	PS1	PS2	PS3	PS4	PS5	PS6
Max Height Range	110µm	300µm	1.1mm	3.5mm	10mm	24mm
Working Distance	3.35mm	10.8mm	12.0mm	16.2mm	25.9mm	20mm
Lateral Accuracy (X-Y)	0.9µm	1.2µm	2.0µm	3.0µm	7.0µm	8.0µm
Height Repeatability (Ra) *	1.2nm	2.2nm	3.4nm	17nm	31nm	41nm

HIGH SPEED SENSOR (192 Points)	LS1	LS2	LS3
Max Height Range	200µm	0.95mm	3.9mm
Working Distance	5.3mm	18.5mm	41mm
Height Repeatability (Ra) *	14nm	21nm	70nm
Line Width	0.96mm	1.91mm	4.78mm
Pitch	5µm	10µm	25µm
Lateral Accuracy of each point	1µm	2µm	5µm
Acquisition Rate (points per second)	384KHz	384KHz	384KHz

* Fixed point on glass, average height variation for 1200points (100 sampling)

Specifications are continuously improving, please contact **NANOVEA** for the latest.

N Today's Standard For Tomorrow's Materials.

NANOVEA instruments can be found in renowned education and industrial [organizations around the world](#).

From [aerospace applications](#) to [medical devices](#), thousands of clients at the frontiers of the most demanding industries, with no room for error, rely on our instruments' [unmatched accuracy](#) and [technical superiority](#).



[Nanovea.com](#)



[Nanovea linkedin](#)



[Nanovea Films](#)



[Nanovea FB](#)



[@nanoveainc](#)

