

***LUBRICATING EYE DROP COMPARISON***  
——— ***USING THE*** ——  
***NANOVEA T50 TRIBOMETER***



Prepared by  
**Stacey Pereira, Jocelyn Esparza, and Pierre Leroux**



## Introduction


Eye drop solutions are used to alleviate symptoms caused by a range of eye problems. For example, they can be used to treat minor eye irritation (e.g. dryness and redness), delay the onset of glaucoma or treat infections. Eye drop solutions sold over-the-counter are mainly used to treat dryness. Their effectiveness in lubricating the eye can be compared and measured with a coefficient of friction test.

## Importance of Testing Eye Drop Solutions

Dry eyes can be caused by a wide range of factors, for example, computer eye strain or being outdoors in extreme weather conditions. Good lubricating eye drops help maintain and supplement the moisture on the outer surface of the eyes. This works to alleviate the discomfort, burning or irritation and redness associated with dry eyes. By measuring the coefficient of friction (COF) of an eye drop solution, its lubricating efficiency and how it compares to other solutions can be determined.

# Measurement Objectives

In this study, the coefficient of friction (COF) of three different lubricating eye drop solutions was measured using the pin-on-disk setup on the Nanovea T50 Tribometer.



**Eye-drop solution sample on Nanovea T50 Tribometer**

# Measurement Parameters

A 6mm diameter spherical pin made of alumina was applied to a glass slide with each eye drop solution acting as the lubricant between the two surfaces. The test parameters used for all experiments are summarized in Table 1 below.

Please note that the alumina pins and glass slides used in this study were chosen due to their availability at the time of testing. A wide range of materials and lubricants can be used to simulate real life applications or tailor the experiment below.

Test Parameters	
Load	1 N
Duration of test	3 min
Rotational rate	200 rpm
Radius of track	2 mm
Revolutions	600
Ball Diameter	6 mm
Ball Material	Alumina
Lubricant	Samples A, B, C, D
Temperature	23°C (room)

Table 1: Test parameters for COF measurements.



Samples of eye-drop solutions

# Results and Discussion

The maximum, minimum, and average coefficient of friction values for the three different eye drop solutions tested are tabulated in Table 2 below. The COF v. Revolutions graphs for each eye drop solution are depicted in Figures 2-4. The COF during each test remained relatively constant for most of the total test duration. Sample A had the lowest average COF indicating it had the best lubrication properties.

Sample	Max COF	Min COF	Avg COF
A	$0.578 \pm 0.205$	$0.18 \pm 0.277$	$0.317 \pm 0.121$
B	$1.05 \pm 0.02$	$0.121 \pm 0.071$	$0.78 \pm 0.024$
C	$1.09 \pm 0.06$	$0.085 \pm 0.069$	$0.84 \pm 0.103$

Table 2: The measure coefficient of friction values for Samples A, B, and C.






Figure 2: Change in COF over number of revolutions for Sample A.

# Results and Discussion



**Figure 3: Change in COF over number of revolutions for Sample B.**



**Figure 4: Change in COF over number of revolutions for Sample C.**



## Conclusion

In this study we showcase the capability of the Nanovea T50 Tribometer in measuring the coefficient of friction of three eye drop solutions. Based on these values, we show that Sample A had a lower coefficient of friction and therefore exhibits better lubrication in comparison to the other two samples.

Nanovea Tribometers offers precise and repeatable wear and friction testing using ISO and ASTM compliant rotative and linear modules. It also provides optional high temperature wear, lubrication, and tribo-corrosion modules available in one pre-integrated system. Such versatility allows users to better simulate the real application environment and improve fundamental understanding of the wear mechanism and tribological characteristics of various materials.

Learn more about [Nanovea Tribometer](#) or [Lab Services](#).

## T2000 Tribometer



### Load Control, High Speed, & High Force

Dual Controlled Loads

Extreme Robust Design

High Speed Sensor Capability

Friction, Wear, Scratch Testing

[Learn More about the T2000](#)

## T50 Tribometer



### Versatile Wear & Friction Tester

Multi-Module System

Speed Control from 0.01-5000 RPM

Robust with Open Platform

Wide Range of Environmental Conditions

[Learn More about the T50](#)

## Thank you for reading!

We appreciate your interest in our technology and services. Read more about all of our product line and lab services at [www.nanovea.com](http://www.nanovea.com)

Call to Schedule a demo today!

If you have any questions please email us at [info@nanovea.com](mailto:info@nanovea.com)

## Recommended Reading


Check out our other application note :

Nano Mechanical Characterization  
of Spring Constants

**NANOVEA**  
A Better Measure

www.nanovea.com  
info@nanovea.com  
euroinfo@nanovea.com  
mexinfo@nanovea.com  
(949)-461-9292

**NANO MECHANICAL CHARACTERIZATION**  
— OF —  
**SPRING CONSTANTS**



Prepared by  
**Xavier Herrera-Keehn & Jocelyn Esparza**

Today's standard for tomorrow's materials. © JULY 2019 NANOVEA