

NANOVEA

DENTAL TOOLS

***DIMENSIONAL & SURFACE ROUGHNESS
ANALYSIS***



Prepared by

Duanjie Li, PhD., Jonathan Thomas, and Pierre Leroux



INTRODUCTION

Having precise dimensions and optimal surface roughness are vital to the functionality of dental screws. Many dental screw dimensions require high precision such as radii, angles, distances, and step heights. Understanding local surface roughness is also highly important for any medical tool or part being inserted inside the human body to minimize sliding friction.

NON-CONTACT PROFILOMETRY FOR DIMENSIONAL STUDY

Nanovea 3D Non-Contact Profilers use a chromatic light-based technology to measure any material surface: transparent, opaque, specular, diffusive, polished or rough. Unlike a touch probe technique, the non-contact technique can measure inside tight areas and will not add any intrinsic errors due to deformation caused by the tip pressing on a softer plastic material. Chromatic light-based technology also offers superior lateral and height accuracies compared to focus variation technology.


Nanovea Profilers can scan large surfaces directly without stitching and profile the length of a part in a few seconds. Nano through macro range surface features and high surface angles can be measured due to the profiler's ability to measure surfaces without any complex algorithms manipulating the results.

MEASUREMENT OBJECTIVE

In this application, the Nanovea ST400 Optical Profiler was used to measure a dental screw along flat and thread features in a single measurement. The surface roughness was calculated from the flat area, and various dimensions of the threaded features were determined.



Sample of dental screw analyzed by NANOVEA Optical Profiler.






Dental screw sample analyzed.

RESULTS

3D Surface

The 3D View and False Color View of the dental screw shows a flat area with threading starting on either side. It provides users a straightforward tool to directly observe the morphology of the screw from different angles. The flat area was extracted from the full scan to measure its surface roughness.






ISO 25178


Height Parameters

Sa	0.9637	μm
Sq	1.220	μm
Sp	5.644	μm
Sv	6.051	μm
Sz	11.70	μm
Ssk	0.2352	
Sku	3.218	


2D Surface Analysis

Line profiles can also be extracted from the surface to show a cross-sectional view of the screw. The Contour Analysis and step height studies were used to measure precise dimensions at a certain location on the screw.





Parameters	Step 1	Unit
Maximum height	0.5724	mm
Mean height	0.5611	mm





CONCLUSION

In this application, we have showcase the Nanovea 3D Non-Contact Profiler's ability to precisely calculate local surface roughness and measure large dimensional features in a single scan.

The data shows a local surface roughness of $0.9637\ \mu\text{m}$. The radius of the screw between threads was found to be $1.729\ \text{mm}$, and the threads had an average height of $0.413\ \text{mm}$. The average angle between the threads was determined to be 61.3° .

The data shown here represents only a portion of the calculations available in the analysis software.