


Quality Analysis on Electrical Discharge Machined Metals

Electrical discharge machining, or EDM, is a manufacturing process that removes material via electrical discharges [1]. This machining process is generally used with conductive metals that would be difficult to machine with conventional methods.

As with all machining processes, precision and accuracy must be high in order to meet acceptable tolerance levels. In this application note, the quality of the machined metals will be assessed with a Nanovea 3D non-contact profilometer.



Importance of Profiling Machined Metals


As machining becomes more precise, simple measurement tools, such as rulers or calipers, no longer meet the level of accuracy needed to satisfy tighter tolerance requirements. High resolution measurement tools must be used. In order to not bottleneck the production process, the measurement time of quality control instruments must match the manufacturing speed as well. The ideal quality control instrument will have high scanning resolution, short measurement time (including sample preparation and setup) and simple user-friendly analysis tools and software.

The Nanovea 3D line sensor is the ideal instrument for quality control inspections of machined metals. It is able to quickly scan and produce high resolution surface scans within seconds which can then be analyzed to quantify multiple dimensions of a machined part(s). Its ability to measure all types of materials, regardless of curvature, makes it ideal for use in measuring a vast range of surfaces in quality control applications.

Measurement Objectives

Equipment Featured

NANOVEA HS2000



High Speed Inspection & Precision Flatness Measure

Advanced Automation

Customizable Options

High Speed


Precision Flatness Measurement

Rigid and Stable Structure

[Learn More about our Profilometers!](#)

Measurement Objectives

Two metallic samples with hexagonal patterned holes drilled via wire EDM were inspected with Nanovea's high-speed line sensor. The roundness, equivalent diameter, and spacing were quantified to observe the precision and accuracy of the WEDM machining process.




Measurement Parameters

Table 1:

Test Parameter	Value
Instrument	HS 2000L
Optical Sensor	LS2 (1000µm height range)
Scan size (mm)	40mm x 20mm
Step size (µm)	10µm x 10µm
Scan time (h:m:s)	00:00:38

Samples Tested



Two metallic samples with hexagonal patterned holes drilled via wire EDM

Profilometry Results

Results: Sample A



Sample A




Figure 1: 3D view of scan taken for Sample A






Figure 2: A) Motif analysis of holes B) histogram of equivalent diameter distribution C) histogram of roundness distribution for Sample A

Profilometry Results



Sample A

Figure 3: Extracted area of hole with outlying roundness for Sample A





Figure 4: Contour analysis of spacing between holes for Sample A

Profilometry Results

Results: Sample B



Sample B

Figure 5: 3D view of scan taken for Sample B





Figure 6: A) Motif analysis of holes B) histogram of equivalent diameter distribution C) histogram of roundness distribution for Sample B

Profilometry Results



Sample B

Figure 7 :Extracted area on outlier in equivalent diameter for Sample B




Figure 8: Contour analysis of spacing between holes for Sample B


Profilometry Results

Discussion

A 40mm x 20mm area was scanned, with nanometer level resolution, within 38 seconds. The results for Sample A and B are very similar. A hole density of approximately 4.25 holes/cm² was seen on both samples. Sample specific details are stated below.

The holes on Sample A were observed to have an equivalent diameter of 0.9887 ± 0.0048 mm and roundness of 0.9630 ± 0.0176 (note: a roundness value of 1 is considered a perfect circle) using a Motifs study. The histograms show that the majority of holes had an equivalent diameter near the mean of 0.9887mm. Two holes are seen to be two standard deviations above the mean. Depending on the acceptable tolerances, these holes may not pass quality control requirements. For roundness, only one hole is significantly different. The outlier can be closely inspected by extracting it and its surrounding area as seen in Figure 3. The average spacing between each hole and its nearest neighbors (3 for hexagonal patterns) is 4.107 ± 0.013 mm.

A Motifs analysis on Sample B shows that the holes on that sample have an equivalent diameter of 0.9842 ± 0.0057 mm and roundness of 0.9630 ± 0.0079 . One hole was singled out from the rest due to its outlying equivalent diameter, which can be seen in Figure 7. The contour analysis in Figure 8 shows the spacing between holes to be 4.107 ± 0.16 mm, almost the same as Sample A.



Conclusion

The equivalent diameter, roundness, and spacing of samples machined via wire EDM were quantified to observe the precision and accuracy of wire EDM machining. Very small differences were found between both samples. Few holes from the two samples measured fell outside several standard deviation from the mean, and were singularly inspected by extracting the area that contained the outlying hole of interest from the original scan. Multiple dimensions were measured on each sample with our analysis software, a software that is capable of creating templates for quick and automatic generation of results like those seen in this document. Our study has shown that Nanovea's profilometry technology is capable of providing quick and accurate measurements for quality control applications on machined parts.

You can explore our extensive digital library of application notes for more applications that our Profilometers have been used for.

<https://nanovea.com/app-notes/profilometry>

References

[1] Jameson, Elman C. Electrical discharge machining. Society of Manufacturing Engineers, 2001.

Thank you for reading!

We appreciate your interest in our technology and services. Read more about all of our product line and lab services at www.nanovea.com

[Contact us today!](#)

If you have any questions please email us at info@nanovea.com

Recommended Reading

Check out our other application note where we conduct a Viscoelastic Analysis on Rubber with Nanoindentation

[Click here to read!](#)

NANOVEA

A Better Measure

www.nanovea.com
info@nanovea.com
euroinfo@nanovea.com
mexinfo@nanovea.com
+1 (949) 461-9292



Viscoelastic Analysis of Rubber with Nanoindentation DMA

Viscoelasticity is referred to as the property of materials that exhibit both viscous and elastic characteristics when undergoing deformation.

A viscous material resists shear flow and strains linearly with time when a stress is applied, unlike an elastic material that strains immediately when stressed and returns to original state once the stress is removed. A viscoelastic material exhibits elements of both properties and therefore has a complex modulus.